

STORAGE OF REAGENTS, CONTROLS AND SUPPLIES

Reagent	UNOPENED	Opened Room Temp	Opened Refrigerator Temp	Notes
HIV ANTIBODIES				
Orasure Oraquick H&L Controls	Store: 2-8 °C (36°-46 °F) Expiration date: Listed	DO NOT STORE AT ROOM TEMPERATURE	Expires in 21 days when stored (2-8 °C) (36°-46 °F)	Do not freeze
Orasure Oraquick HIV-1 <i>Classic</i>	Listed Expiration Date when stored at: 2-27 °C (35°-80 °F)	Listed Expiration Date when stored at: 2-27°C (35-80°F) no sunlight	Do not freeze Do not store opened!!	BE CAREFUL!! <ul style="list-style-type: none"> • DO NOT STORE >80° F • RUN QC IF STORAGE WAS OUT OF RANGE
Orasure Oraquick HIV-1/2 <i>Advanced</i>	Listed Expiration Date when stored at: 2-27 °C (35°-80 °F)	Listed Expiration Date when stored at: 2-27°C (35-80°F) no sunlight	Do not freeze Do not store opened!!	BE CAREFUL!! <ul style="list-style-type: none"> • DO NOT STORE >80° F • RUN QC IF STORAGE WAS OUT OF RANGE